

# Liviano Carbon Fiber Rollator

ULTRA-LIGHTWEIGHT CARBON FIBER DESIGN



## Product Specifications

| Model              | Liviano Rollator |
|--------------------|------------------|
| Material           | Carbon Fibre     |
| Weight Capacity    | 330 lbs (150 kg) |
| Max Carry Bag Load | 13 lbs (6 kg)    |
| Total Weight       | 12 lbs (5.5 kg)  |
| Seat Width         | 18"              |

|               |                                           |
|---------------|-------------------------------------------|
| Handle Height | 27" - 33.5" (Standard) 31" - 36.5" (Tall) |
| Seat Height   | 22" (Standard) 25" (Tall)                 |
| Wheel Size    | 8"                                        |
| Unfolded Size | 25.5" (D) x 24.5" (W) x 30.7" (H)         |
| Folded Size   | 25.5" (D) x 8.6" (W) x 32.5" (H)          |
| Package Size  | 26" (D) x 9" (W) x 34.6" (H)              |

## Key Features



### Ultra-Lightweight:

At just 12 lbs, the Liviano is one of the lightest rollators available — ideal for everyday use, car travel, or even taking on a plane.



### Compact Folding:

Folds down to only 8.5" wide, making it easy to store in tight spaces like car trunks, closets, or under the bed.



### One-Handed Fold:

A smooth upward lift of the seat handle folds the rollator quickly and easily.



### Strong & Stylish:

Crafted from durable carbon fiber for maximum strength without bulk, while offering a sleek modern look.



### Carry Basket:

Convenient storage basket supports up to 13 lbs, perfect for groceries, essentials, or personal belongings.




### Comfortable Backstrap:

Built-in backstrap provides reliable posture support whenever you stop to rest.

## Colour Options:



Height Adjustable  
Handles



Latch keeps  
rollator closed



Extra Large  
wheels



Carbon-fiber  
lightweight frame



Ruby  
Red

Onyx  
Black

