

SolCare 200 Mattress

PREMIUM PRESSURE REDISTRIBUTION MATTRESS

WWW.TEKNOCAREMEDICAL.CA



Product Specifications

Mattress Model	SolCare 200
Dimensions	80" (L) × 36" (W) × 6" (H)
Weight Capacity	350 lbs (158 kg)
Mattress Cover Material	4-way stretch nylon
Warranty	4 years foam 2 years cover

Key Features



Multi-Tier Foam Structure:

Even pressure redistribution for high-risk areas.



High-Density Foam Layer:

Enhanced comfort for torso & head.



Memory Foam Heel Zone:

Superior cushioning, reduces bottoming out.



Air-Channel Foam with temperature control:

Combines Vapor-permeable fabric and open-cell foam to improve airflow and maintain a comfortable sleeping temperature



Reinforced Edges:

Safer transfers and bedside sitting.



4-Way Stretch Cover:

Water-resistant, antimicrobial, low-shear.



Protective Zipper Flap:

Guards against fluids & contamination.



Durable Build:

Long-term durability for at home or facility care

