

VitalFlex Eco

AFFORDABLE, FUNCTIONAL
& SAFE HOME HOSPITAL BED

WWW.TEKNOCAREMEDICAL.CA



Product Specifications

Bed Model	VitalFlex Eco
Product Size	36" x 84"
Sleep Deck Size	36" x 80"
Product Weight	80 lbs
Load Capacity	450 lbs
Height Range	9"-24"
Max Backrest Angle	65°
Max Legrest Angle	35°
Cardiac Chair Position	Yes



Backrest
0°-65°



Footrest
0°-35°



Height
adjustment
9" - 24"

Key Features



Ultra-Low Safety:

Lowers to 9" to help prevent fall-related injuries.



Cardiac Chair & Tilt:

One-touch cardiac chair and reverse Trendelenburg for therapeutic comfort.



Adjustable Comfort:

Backrest up to 65° and footrest up to 35°.



Durable Deck:

Metal grid mattress deck for stability and airflow.



Night Safety:

Built-in underbed lighting for nighttime use.



Convenience:

Dual USB chargers for devices and ergonomic 10-position handset.



Emergency Backup:

9V battery backup included.



Optional Accessories:

Half siderails or rotating assist handle available.

Siderail Options



Half Rails



Rotating Assist Handle

