

VitalFlex Elite Care Bed

COMFORT, SAFETY & STYLE IN ONE BED


WWW.TEKNOCAREMEDICAL.CA



Backrest
0°-70°



Footrest
0°-30°



Front & Reverse
Trendelenburg
0-12°




Height
9" - 31"

Product Specifications

Bed Model	VitalFlex Elite
Sleep Deck Dimensions	36" x 80"
Overall Bed Dimensions	41" x 84"
Height Range	9" - 31"
Patient Weight Capacity	420 lbs
Battery Backup	Included

Key Features



Height Adjustment:

Electrically adjustable from 9" (ultra-low) to 31" (care height) for maximum safety and flexibility.



Backrest & Legrest:

Backrest tilts up to 70°, legrest adjustable to 35° for patient comfort and caregiver ease.



Comfort Positioning:

Pre-programmed "cardiac chair" setting with one-touch ergonomic handset (locking key for safety).



Design & Finish:

Solid wood paneling in two tones for a modern, homelike look that blends with interiors.



Mobility & Safety:

Central parking brake, smooth maneuverability with 4x 3" casters + 2x 4" casters.



Cleaning & Maintenance:

Durable powder-coated metal components; easy to clean, long-lasting build.

Optional Side Rails



Folding Steel
Full Rail



Steel & Wood
Half Rail

Board Colours



Walnut

Oak

Accessories



Extension
to 84"



Underbed
Lighting

