

SolCare 100 Mattress

PREMIUM THERAPEUTIC MEDICAL MATTRESS

WWW.TEKNOCAREMEDICAL.CA



Product Specifications

Mattress Model	SolCare 100
Dimensions	80" (L) × 36" (W) × 6" (H)
Weight Capacity	350 lbs (158 kg)
Mattress Cover Material	4-way stretch nylon
Warranty	4 years foam 2 years cover

Key Features



Designed for Low to Moderate Risk Patients:

Provides dependable comfort and pressure redistribution.



Ultra-Soft Foam Formulation:

Allows deep immersion, reducing pressure on sensitive areas.




Low-Shear, Antimicrobial & Fluid-Resistant Material:

Protects skin integrity and minimizes friction.



Temperature Control:

Combines Vapor-permeable fabric and open-cell foam to improve airflow and maintain a comfortable sleeping temperature



High-Resilient Foam Core:

1.5 lb density foam ensures durability and long-lasting support.



Protective Nylon Cover:

Vapor-permeable top cover minimizes perspiration and keeps surface cool and dry.



360° Zipper Design:

Easy removal for cleaning or replacement.



High density Foam Perimeter:

For edge support and easier transfers

