

Vital Lite Wheelchair

LIGHTWEIGHT, PORTABLE & DURABLE

WWW.TEKNOCAREMEDICAL.CA



Product Specifications - VL-WC-16 / VL-WC-18

Model	Vital lite Wheelchair
Weight Capacity	220 lbs.
Overall Weight	30 lbs.
Seat Width	16" , 18"
Seat Depth	16"
Seat-to-Floor Height	19"
Overall Dimensions Open	36" L x 25"/27" W x 35" H

Overall Dimensions Closed	36" (L) x 12.7" (W)
Footrests	Swing-Away
Rear Wheel Size	22"
Front Caster Size	7"
Carton Shipping Weight	38 lbs
Shipping Dimensions:	30.7" x 13.4" x 32.2"

Key Features



Lightweight & Portable:

Weighs only 30 lbs, making it easy to handle, store, and transport.



Compact Folding Design:

Backrest folds down and frame collapses for easy storage in tight spaces or car trunks.



Flip-Up Armrests:

Full-length, padded armrests that flip back for easier transfers and added convenience.



Strong Frame:

Built with a rust-resistant, powder-coated aluminum alloy frame for durability and long-lasting use.



Swing-Away Footrests:

Detachable, swing-away footrests improve accessibility and support.



Comfortable Seating:

Cushion made from breathable Oxford fabric with mesh for all-day comfort.



Safe & Reliable:

Equipped with self-locking nursing brakes, wheel locks, and a double-cross axle for enhanced safety and stability.



User-Friendly Rims:

Ergonomic hand rims are designed for smoother propulsion and less strain on hands.

Colour Options:



Flip Back
Armrests



Dual Brake
System



Quick-release
footrests



Ergonomic
22" Wheels



Seafoam
Green Midnight
Black

