

Vital Lite Transport Chair

COMPACT, LIGHTWEIGHT & CARE-FRIENDLY

WWW.TEKNOCAREMEDICAL.CA



Product Specifications - VL-TC-16 / VL-TC-18

Model	Transport Chair
Weight Capacity	220 lbs.
Overall Weight	21 lbs.
Seat Width	16" or 18"
Seat Depth	16"
Seat-to-Floor Height	19"
Dimensions Open with Footrest:	34" L x 22.5"/24.5" W x 34.5"H

Dimensions Closed, Footrests Up:	26" (L) x 12" (W)
Footrests	Flip Back
Rear Wheel Size	15" MAG PU tire
Front Caster Size	7"
Carton Shipping Weight	26 lbs
Shipping Dimensions:	28.7" x 11.8" x 31.5"

Key Features



Ultra-Lightweight:

Weighs only 21 lbs, making it easy for caregivers to maneuver and transport.



Compact Folding Design:

Backrest and frame fold down for quick storage in tight spaces or car trunks.



Flip-Up Armrests:

Full-length padded armrests flip up for easier side transfers.



Durable Frame:

Constructed with a rust-resistant, powder-coated aluminum alloy frame for strength and reliability.



Innovative Footrests:

Detachable, flip-back adjustable footrests enhance comfort and accessibility.



Comfortable Seating:

Breathable Oxford and mesh fabric cushion ensures comfort during extended use.



Safe & Stable:

Equipped with nursing brakes, wheel locks, and a double-cross axle for safety and smooth handling.



Effortless Maneuverability:

15" rear MAG PU tires and 7" front casters ensure smooth indoor and outdoor navigation.

Colour Options:



Flip Back
Armrests



Larger 16"
Rear Casters



Flip Back
Footrests



Dual Brake
System



Seafoam
Green

Midnight
Black

