

VitalFlex Standard & VitalFlex Low

FLEXIBLE, RELIABLE & SAFE CARE BEDS



Backrest
0°-70°



Footrest
0°-30°



Height
Std: 14"-32"
Low: 8"-25"



Front & Reverse
Trendelenburg
0°-12°

Product Specifications

Bed Model	VitalFlex Standard	VitalFlex Low
Sleep Deck Dimensions	36" x 80"	36" x 80"
Overall Bed Dimensions	40" x 85"	40" x 85"
Height Range	14" - 32"	8" - 25"
Total Bed Weight	200 lbs	200 lbs



Key Features



Height Adjustment:

VitalFlex Standard adjusts from 14"-32"; VitalFlex Low lowers to 8"-25" for maximum safety.



Bed Platform:

4-section steel slat deck (36" x 80") with integrated mattress guides at head & foot.



Power & Emergency:

24V power supply for greater safety and 9V battery backup for emergencies.



Design Options:

VitalFlex available with Oak or Walnut headboards.



Backrest & Footrest:

Electrically adjustable backrest (0°-70°) and footrest (0°-30°) for patient comfort.



Safety & Mobility:

Mechanical foot elevation (5 positions), 1.5" x 4" individually lockable casters.



Optional Accessories:

Bed extension up to 8", patient lifting pole (supports up to 176 lbs), optional siderails (foldable steel, split with wood top, assist rail).



Compact Packaging:

Easy-to-handle packaging for shipping and setup.

Optional Side Rails



bed Assist Rail
Gab Handle



Solid Wood Top
Half Rails

Wood Headboard Colour



Oak



Walnut

