

Symphony Hi-Low Bed

WWW.TEKNOCAREMEDICAL.CA



Backrest
0°-58°



Footrest
0°-38°



Front & Reverse
Trendelenburg
0-12°



Height
adjustment
10.5" - 24"



German
Brand Okin
motor

Product Specifications

Bed Model	Twin XL	Queen	Split King
Product Size	38" x 80"	60" x 80"	76" x 80"
Product Weight	198 lbs	231 lbs	396 lbs
Weight Capacity	400 lbs	400 lbs	800 lbs
Height Range	10.5"-24" (no legs or mattress)		

Package Size

Bed Model	Carton Dimension	Gross Weight
Twin XL	82" x 41" x 11"	198 lbs
Queen	2x of 42" x 63" x 8"	143 lbs & 88 lbs
Split King	2x 82" x 41" x 11"	2x 198 lbs

Bed Leg Option

Plastic Leg Height	3.5" (2.4" diameter)
2" Casters Height	2.4" up to 800 lbs weight

Key Features



Height Adjustment:

Electrically adjustable from 9"-24" for versatile use and ease of care.



Trendelenburg Function:

Trendelenburg & Reverse Trendelenburg up to 12° for therapeutic support.



Massage Function:

Two-zone massage for comfort and relaxation.



Wireless Remote:

Intuitive handset with multiple options to find ideal comfort.



Fits Your Style:

Includes bracket that fits any traditional headboard.



Backrest & Footrest:

Backrest tilt up to 58°, footrest tilt up to 38° for optimal positioning.



Underbed Lighting:

Built-in safety lighting for nighttime use.



Smart Connectivity:

Integrated Bluetooth and USB outlets for device charging.



Motor:

German Brand Okin motor for reliable, smooth performance.



A Great Night's Rest:

Packaged mattress features 2" gel-infused memory foam and -6" high-density foam core

Accessories



Voice Activated controls with wireless phone charger



Height Adjustable Side Rails

